

Destiny™

[with BluWave™ SP]



Destiny:

Congratulations
on your journey
to improved
hearing.

Only **one** out of **five** people with hearing loss actually **seeks help**, so by visiting your hearing care professional, you've already taken a big step. And with **Destiny**, you'll be wearing one of the **world's most advanced hearing aids**: fully automatic, comfortable, and crystal clear.

The smartest hearing science

Only Destiny has BluWave Signal Processing, Starkey's proprietary system designed to optimize noise management, feedback and whistling cancellation and talking on the phone with no trouble.

The Destiny series includes multiple technology levels. Each comes with our best-in-class Active Feedback Intercept, which virtually eliminates the nuisance of feedback. Your hearing care professional will help determine the levels of environmental adaptation and directional capabilities you need.

The Destiny 1600 is our most automatic hearing instrument on the market today, featuring touchless operation and a host of advanced technology benefits, including Self Check, Voice Indicators, and Follow-Up Reminder.



BTE



OTE



mini



**HERE'S WHAT HEARING AID WEARERS
ARE SAYING ABOUT DESTINY:**

"This is a miracle."

"I'm understanding every word."

"No feedback at all — that's amazing."

"I can't imagine anything better than Destiny."



ITE



ITC



CIC

SMART FEATURES

- Touchless
- Voice Indicators
- Self Check
- Reminder



NOISY ENVIRONMENTS


- In the stands
- At a party
- On the town
- Out of doors



- At worship



- In a meeting



- A quiet café




- One-on-one



- At home



QUIET ENVIRONMENTS



Hear in style

You'll find there are a variety of hearing aid styles from which to choose. Technological advances have resulted in more options to meet specific hearing problems and comfort preference. The choice of technology in the hearing aid is primarily what makes the difference in performance and price, but style also plays a role.

THERE ARE SIX HEARING AID STYLES.



Behind the Ear (BTE)

- *Easiest to adjust*
- *Recommended for most types of hearing loss and ages*
- *Attaches to earbud or custom earmold*



mini

- *Smallest BTE model*
- *Available with thin tubing and tiny earbud for ultimate discretion*
- *Appropriate for most types of hearing loss and ages*



On the Ear (OTE)

- *Lightweight with a micro-thin tube and tiny earbud for ultimate discretion*
- *Appropriate for mild to moderately-high frequency hearing loss*



In the Ear (ITE)

- *Easy to adjust*
- *Appropriate for mild to severe hearing loss*
- *Custom-made case for comfort*



In the Canal (ITC)

- *Less visible*
- *Appropriate for mild to mildly-severe hearing loss*
- *Variety of circuit choices that previously required a larger aid*



Completely in the Canal (CIC)

- *Virtually invisible*
- *Appropriate for mild to moderately-severe hearing loss*

Destiny Features

Active Feedback Intercept (AFI)*

The most effective feedback cancellation system on the market today. AFI uses an entirely new method that effectively erases whistling, squealing, and other irritating feedback sounds.

Directional Speech Detector (DSD)

Standard on behind-the-ear styles and optional on many others, DSD offers superb directionality by focusing on speech sounds and reducing background noise. This technology continuously monitors environments and adjusts accordingly.

Acoustic Signature

A state-of-the-art sound management tool available on the 1200 and 1600 series. This system detects and classifies different sound environments, then adapts automatically to an optimal setting. Move from crowded restaurant to taxi to home, without ever making a manual adjustment.

Smart Style

Our commitment to making Destiny not only comfortable, clear and supremely effective, but attractive as well. Our Destiny BTE, mini and OTE come in a range of sleek colors designed to match your hair color, skin tone, or lifestyle. And our range of styles and sizes means you can find an option that looks and feels natural on you.

DESTINY ALSO OFFERS A FULL PACKAGE OF SMART FEATURES, INCLUDING:

Voice Indicators

The option of using speech, with male or female voices, to alert you to low battery, memory, as well as others.

Self Check

Simply opening and closing the battery door three times runs a diagnostic check of the microphone, circuit and receiver.

Reminder

To help remind you of maintenance and return visits, we've included an audible follow-up reminder.

* Independent studies prove that AFI is the very best in its class. Our patent-pending method for feedback cancellation uses a sophisticated sub-band system with 16 pathways for incoming sound and provides better gain margin improvement than any other hearing instrument on the market today.



A global hearing technology company
headquartered in Eden Prairie, Minnesota, U.S.A.

Starkey Laboratories, Inc.
6700 Washington Avenue South
Eden Prairie, MN 55344-3476
800.328.8602

www.starkey.com

© 2006-2008 Starkey Laboratories, Inc. All Rights Reserved
79522-000 03/08 BROC9522-06-EE-ST